

COMPACT PRECISION SCALE UC-322



Precision Accuracy in a Compact Size 8" x 11"

The UC-322 personal weight scale is an accurate method of tracking your weight loss and body mass index. With 0.1 lb increments, it is easy to notice even the smallest weight loss. Sized smaller than an average sheet of paper, the UC-322 has a capacity of 330 lb. to accommodate people of all sizes and weights.

The UC-322 provides dual memory zones with 15 memory recall each. By using the total 30 memory capacity, two individuals can track over 2 weeks of weight trend (one measurement per day).

Automatically calculates Body Mass Index (BMI) according to weight-to-height ratio.

BMI	Weight Status
Below 18.5	Underweight
18.5 - 24.9	Normal
25 - 29.9	Overweight
30 and above	Obese

- Accurate readings within 0.1 LBS.
- Compact size with 330 LB. capacity
- Body mass index (BMI) calculation
- Dual Memory Zones with 15 memory recall each
- Low profile for ease of use



For questions about LifeSource products, call us toll-free at 1-800-726-3364 and a specially trained representative will assist you, or visit us at www.LifeSourceOnline.com.

UC-322 Compact Precision Scale

PART NUMBER	UPC NUMBER	Capacity x Resolution
UC-322	0-93764-60202-3	330 lb x 0.1 lb

This compact precision scale includes:

- Instruction manual
- Four "AA" batteries
- Four elevation feet

Technical Specifications

Capacity	330 lb (150 kg)
Resolution	0.1 lb (0.05 kg)
Dimensions	11.5" (w) x 8" (d) x 1.25" (h) (295 x 210 x 40 mm)
Weight	Approx. 3.5 lb (1.6 kg)
Display	LCD, character height: 1 inch (25 mm)
Sensor	Four Load Cells
Power Source	Four type "AA" batteries included
Battery Life	Approx. 2,000 measurements
Auto Power Off	5 seconds
Operating Temperature	50 - 95° F (10 - 35° C) RH 30% - 85%

Manufactured by an ISO 9001 and ISO 13485 registered company. Engineered in Japan. Made in China.

Packaging and Shipping Specifications

Selling Unit (Box)	Dimensions (L x W x H) 8.25" x 13.75" x 3.38"	Weight 5 lb	Individual Item Cube 0.22 Cubic Feet
Unit Case (5 Units)	Dimensions (L x W x H) 14.5" x 18" x 12"	Weight 25 lb	Case Cube 1.8 Cubic Feet



1756 Automation Parkway
San Jose, CA 95131
1-888-726-3364
www.LifeSourceOnline.com



Ask your A&D Sales Representative about other LifeSource products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Stethoscopes

For questions about LifeSource products, call us toll-free at 1-800-726-3364 or visit our web site at www.LifeSourceOnline.com